

A Novel Method For Discovery of Peripheral Blood Biomarkers in Idiopathic Pulmonary Fibrosis Using Extensive Depletion and TMTcalibrator™ Tissue-Enhanced Plasma Proteomics

I. Pike¹, M. Bremang¹, P.J. Wolters², R. Gasier³, S. Turner³, S. Turner³, M. Decaris³
 1 – Proteome Sciences plc, Hamilton House, Mableton Place, London, WC1H 9BB, UK; 2 – University of California at San Francisco, San Francisco, CA 94063, USA;
 3 – PLIANT Therapeutics, 700 Saginaw Drive, Suite 150, Redwood City, CA 94063, USA

Introduction
 Peripheral biomarkers related to the pathogenesis of idiopathic pulmonary fibrosis (IPF) are urgently needed to improve diagnosis, selection and assessment of treatment, particularly in the context of new drug development. Recent improvements in the depletion of high and medium abundant proteins and development of tissue-enhanced fluid proteomics have increased the breadth and depth of plasma proteome coverage. We have now combined these methods to obtain unparalleled coverage of the IPF plasma proteome.

Method
 Longitudinal archival plasma samples (n=25) from six individuals with IPF and five control individuals were obtained from the UCSF Biobank and processed as shown in Figure 1.

- Super-depletion of the top ~70 high and medium abundant proteins using Seppro® IgY14 and Supermix® columns (Merck).
- Plex #1-5 each comprised 6 plasma samples with 4 tissue trigger/calibrant samples. Plex 6 comprised 10 plasma samples only.
- Each TMT® 10plex was analysed using 30 x 2 hour gradient on an Orbitrap® Fusion Tribrid® mass spectrometer with in-line uHPLC.
- MS files were processed with a sequential SEQUEST search strategy in Proteome Discoverer v1.4 (Thermo Scientific) with a standard search followed by a bespoke method in the residual unmatched spectra with variable modification for less-common post-translational modifications such as hydroxylation of proline and lysine.
- Data integration, normalisation, quantitation and statistical analyses were performed within Proteome Sciences' in-house bioinformatics workflows.

• Reported plasma abundance of select IPF-related proteins was obtained through the human plasma proteome database and/or PubMed searches.




Figure 1 – TMTcalibrator™ workflow with plasma Super-depletion and incorporation of four-point trigger/calibrant. TMT126 – 131N = Tandem Mass Tag reporter ion mass-to-charge.

Gene Symbol	TMTcal+SD (# peptides)	SD (# peptides)	Plasma Conc'n
POSTN	195	162	200 ng/ml
Collagen IIV	1240	371	~100 ng/ml
CHI3L1	10	11	40 ng/ml
CCL18	3	0	31 ng/ml
SFA	70	31	30 ng/ml
TGFBI	8	20	5 ng/ml
KRT18	24	0	4 ng/ml
MMF7	2	1	4 ng/ml
KRT8	42	1	3 ng/ml
TGF2	2	0	0.2 ng/ml
LOXL2	1	0	0.1 ng/ml
TGM2	30	1	NR
ITGAV	55	6	NR
ITGB6	22	0	NR
MUC1	19	2	NR
TGFBF1	1	2	NR

Table 2 – Peptide coverage of select IPF-related proteins. TMTcal+SD = Super-depleted plasma with four-channel trigger/calibrant; SD = Super-depletion only. Concentrations from pubmed/proteinatlas.org; NR = not report ed




Figure 2 – SDS-PAGE gels showing consistent plasma depletion – Left panel = IgY14 depletion; Right panel = Supermix depletion




Figure 3 – Biological coverage of TMTcalibrator™ vs Super-depletion alone. A – GO Cellular Component; B – GO Biological Process. TMTcal+SD = Super-depleted plasma with four-channel trigger/calibrant; SD = Super-depleted plasma only

Unique Gene Names	TMTcalibrator™									
	TMT1	TMT2	TMT3	TMT4	TMT5	TMT6	TMT7	TMT8	TMT9	TMT10
Stops Number	164,217	164,242	154,841	176,938	172,692	181,225				
PSMs in Plasma channel	51,091	100,397	77,479	125,373	118,110	98,220				
PSMs in IgY channel	97,890	76,698	66,954	84,543	81,953	40,938				
Peptides	8,903	10,640	9,735	11,927	10,909	4,617				
Protein Groups	7,535	9,846	8,291	9,228	9,146	3,880				

Table 1 – Output metrics for TMTcalibrator™ and TMTMS2 analysis of IPF plasma – Standard Search (upper panel); Sequential FTM search (lower panel)




Figure 4 – Increased coverage of the proteome using TMTcalibrator™ with Super-depletion. A – Protein class; B – Lung proteins (based on 188 genes defined as overexpressed in lung at proteomics.org). TMTcal+SD = Super-depleted plasma with four-channel trigger/calibrant; SD = Super-depleted plasma

Results & Conclusion

- Super-depletion removed ~99% of total plasma protein (Figure 2).
- Detection limit in mid-pg/ml range in plasma channels when using the tissue trigger.
- ~4,000 proteins quantified with Super-depletion only, more than doubled to ~9,000 when tissue trigger added (Table 1).
- Key post-translationally modified proteins only detected when using tissue trigger (Table 1).
- Improved coverage of select IPF-related proteins (Table 2).

- Distributions for GO cellular components (Figure 3A) or biological processes (Figure 3B) were similar with or without tissue trigger.
- Proportionally greater impact on detection of cellular vs secreted protein classes (Figure 4A) and for lung-associated proteins with TMTcalibrator™ than in Super-depleted plasma alone (Figure 4B).
- TMTcalibrator™ with Super-depletion provides outstanding proteome coverage for plasma biomarker discovery offering deeper insights into IPF disease processes and response to treatment.
- Larger studies required for replication/validation of novel candidate biomarkers.

IP & MB are paid employees and hold stock and/or stock options in Proteome Sciences plc. IP is an inventor on patents covering the TMTcalibrator™ technology. ST & MD are paid employees and hold stock and/or stock options in PLiant Therapeutics Inc, who funded the study.

